## Instructions for Use

BISCO, Inc.  
1100 W. Irving Park Rd.  
Schaumburg, IL 60193  
U.S.A.

1-847-534-6000  
1-800-247-3368

IN-154R12  
Rev. 9/19

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed Dentist.

**BisCem®**

Self-Adhesive Resin Cement

GENERAL INFORMATION

**BisCem** is a self-curing, self-adhesive, dual-cure resin cement that is exclusively formulated for luting crowns, bridges, inlays, onlays and posts (prefabricated metal and non-metal fiber posts, as well as cast posts).

2. Resin crown, bridges, inlays and onlays

3. Metal (prefabricated or cast) and non-metallic endodontic posts

4. Porcelain crowns, inlays and onlays (includes alumina and zirconia)

NOTE: It is recommended to use **CHOICE™ 2** for veneer cementation.

Warnings:

• Saliva and oral fluid contamination will seriously compromise dentin bonding.

Cautions:

• Product may contain items that are designed for one use. Use, disposal of used or contaminated tips. Do not clean, disinfect or reuse.

• Working and Setting Times will be accelerated in the warm intraoral environment.

• Do not use hydrogen peroxide, EDTA, or eugenol containing materials prior to the adhesive procedure as these materials may inhibit bonding.

• Since ceramics and porcelain vary greatly, it is important to consult with the ceramic/porcelain manufacturer or your laboratory for instructions on proper treatment and proper prep design prior to using BisCem.

• Avoid contact with the skin; unpolymerized (meth)acrylic resins may cause skin sensitization in susceptible persons. In case of contact, wash with soap and water.

• To prevent cross-contamination, the Dual-Syringe static mixing tip should be left in place until the next application. If used intraorally, it should be wiped with disinfectant prior to storage.

• Since ceramics and porcelains vary greatly, it is important to consult with the ceramic/porcelain manufacturer or your laboratory for instructions on proper surface treatment.

• See individual component labels for specific expiration dates.

• Safety data sheet available at [www.bisco.com](http://www.bisco.com).

## HELPFUL HINTS

• BisCem left on a mixing pad may not appear to set. This is due to oxygen inhibition of the polymerization process. Rest assured the cement will polymerize (set) in a normal fashion under the restoration.

## TECHNICAL INFORMATION

## INSTRUCTIONS FOR USE

For: Crowns, Bridges, Inlays and Onlays

## Preparation Procedures

NOTE: BisCem is a self-adhesive cement which does not require etching or priming of the prepared surfaces.

1. Remove the temporary restoration and all temporary residue.

2. Try in the final restoration to insure proper fit.

3. Attach a mixing tip to the dual-syringe by lining up the key and keyway. Then turn the brown locking mechanism clockwise.

4. Pressing the plunger will mix and dispense the BisCem. (See "TECHNICAL INFORMATION" Section for Working and Setting Time information.)

## Cleaning Procedures

NOTE: Since ceramics and porcelains vary greatly, it is important to consult with the ceramic/porcelain manufacturer or your laboratory for instructions on proper surface treatment.

1. Apply cement to the internal surfaces of the restoration. For inlays, it may be easier to apply the cement to the tooth.

2. Seal the restorative and remove any excess cement. HINT: To aid in the removal of excess cement, initially light cure the margins for 2-3 seconds.

3. Light cure for 20-30 seconds, or allow to self cure. (When using self-cure mode, light curing the margins even briefly will expedite the final set.)

For: Posts (pre-fabricated, metal, cast and non-metal/fiber posts)

## PREPARATION PROCEDURES

1. Remove the syringe cap or mix tip from the dual-syringe.

2. Dispense a small amount of material onto a mixing pad to eliminate any voids in each chamber of the dual-syringe.

3. Attach a mixing tip to the dual-syringe by lining up the key and keyway. Then turn the brown locking mechanism clockwise.

4. Pressing the plunger will mix and dispense the BisCem. (See "TECHNICAL INFORMATION" Section for Working and Setting Time information.)

## Cement Mixing Procedures

1. Apply the syringe cap or mix tip from the dual-syringe.

2. Seal the post and remove any excess cement. HINT: To aid in the removal of excess cement, initially light cure the margins for 2-3 seconds.

3. Light a tip to the dual-syringe by lining up the key and keyway. Then turn the brown locking mechanism clockwise.

4. Process with built-up procedure.

**HYGIENE:** Use of commonly available hygienic protective covering to avoid any contamination of the BisCem cement dual-syringes during the treatment is recommended.

**DISPOSAL:** Refer to community provisions relating to waste. In their absence, refer to national or regional provisions relating to waste.

**STORAGE:** Refrigeration is necessary when in use (2°C/36°F - 8°C/46°F). Allow refrigerated BisCem to reach room temperature before use.

**HYGIENE:** Bisco, Inc. recognizes its responsibility to replace products if proven to be defective. Bisco, Inc. does not accept liability for any damage resulting, either direct or consequential, stemming from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

\* BisCem is a registered trademark of Bisco, Inc.

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CHOICE is a trademark of Bisco, Inc.