## Instructions for Use

**BISCO, Inc.**  
 1100 W. Irving Park Road  
 Schaumburg, IL 60193  
 U.S.A.  
 1-800-247-3368  
 1-847-534-6000  
[www.bisco.com](http://www.bisco.com)

**EC REP** BISCO France  
 208, allée de la Coudoulette  
 13680 La Ciotat de Provence  
 France  
 33-4-90-42-92-92

IN-168R7  
 Rev. 1/21

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed Dentist.

### en PRO-V COAT® Provisional Separating Material

#### GENERAL INFORMATION

**PRO-V COAT** is a "first of its kind" separating agent that allows the clinician to utilize the Immediate Dentin Sealing technique without the risk of the provisional bonding to the tooth structure. The proprietary formulation is a water-soluble separating agent that serves to prohibit bonding of the provisional to the adhesive interface, protects the dentin surface from contamination and cleans up effortlessly without any residue.

#### Indications for Use:

1. Immediate Dentin Sealing (IDS) Technique prior to Provisionalization
  - A. Inlay/Onlay IDS
  - B. Crown IDS

#### Warnings:

- **PRO-V COAT** is highly flammable.

#### Precautions:

- Do not refrigerate **PRO-V COAT**. Store at room temperature (20°C/68°F - 25°C/77°F).
- See individual component labels for specific expiration dates.
- Avoid contact with the skin; unpolymerized resins may cause skin sensitization in susceptible persons. In case of contact, wash skin with soap and water.
- Safety data sheet available on request.
- Safety data sheet available at [www.bisco.com](http://www.bisco.com).

#### INSTRUCTIONS FOR USE

##### IMMEDIATE DENTIN SEALING TECHNIQUE PRIOR TO PROVISIONALIZATION

###### A. INLAY/ONLAY: IDS TECHNIQUE

1. Isolate and prep tooth per inlay/onlay preparation guidelines.
2. Apply adhesive according to manufacturer's instructions.
3. If necessary, block proximal wall undercuts with a flowable and light cure.
4. Remove air inhibited layer with an alcohol pellet. Use finishing bur to re-prep and smooth occlusal enamel margins.
5. Remove rubber dam and take final impression.
6. Isolate preparation and apply **PRO-V COAT** to lubricate IDS. Gently air dry (from 8-10cm from the prep) for 10-15 seconds to evaporate the solvent.
7. Build-up temporary filling material according to manufacturer's instructions.

###### B. CROWN: IDS TECHNIQUE

1. Isolate and prep tooth per crown preparation guidelines.
2. Apply adhesive according to manufacturer's instructions.
3. If necessary, block potholes and undercuts with a flowable and light cure.
4. Remove air inhibited layer with an alcohol pellet.
5. Take final impression.
6. Isolate preparation and apply **PRO-V COAT** to lubricate IDS. Gently air dry (from 8-10cm from the prep) for 10-15 seconds to evaporate the solvent.
7. Build temporary, re-apply **PRO-V COAT**. Gently air dry (from 8-10cm from the prep) for 10-15 seconds to evaporate the solvent.
8. Cement with a temporary resin cement.
9. Check occlusion and adjust slightly out of occlusion.

###### PROVISIONAL REMOVAL

1. Insert a suitable sharp instrument (probe/scaler) into the material and remove in the line of draw.
  2. Thoroughly clean the preparation (with pumice slurry or air abrasion).
  3. Proceed with permanent bonding protocol.
- DISPOSAL:** Refer to community provisions relating to waste. In their absence, refer to national or regional provisions relating to waste.
- STORAGE:** Store **PRO-V COAT** at room temperature (20°C/68°F - 25°C/77°F). Do not refrigerate **PRO-V COAT**.

**WARRANTY:** BISCO, Inc. recognizes its responsibility to replace products if proven to be defective. BISCO, Inc. does not accept liability for any damages or loss, either direct or consequential, stemming from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

\* **PRO-V COAT** is a registered trademark of BISCO, Inc.

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