

Bisco

CE 0459

BISFIL 2B *Self-Cured*

Base Increment Composite

Instructions for Use



BISCO

BISCO, Inc.
1100 W. Irving Park Rd.
Schaumburg, IL 60193
U.S.A.

1-847-534-6000
1-800-247-3368

IN-0041R11
Rev. 11/17

BISFIL™* 2B

Base Increment Composite

GENERAL INFORMATION

BISFIL 2B is a self-cured, radiopaque composite intended for use as the base increment of posterior composite restorations in the Directed Shrinkage Technique. This technique suggests that in contrast to light-cured composites, **BISFIL 2B** begins curing adjacent to the naturally warmer surfaces of the dentin. Since all composites shrink upon polymerization, the net effect of the Directed Shrinkage Technique is to direct polymerization toward the tooth, improving adaptation to the preparation.

BISFIL 2B is a low viscosity, syringeable composite available in shades A3/A3.5 and Universal (original). It is intended for use with adhesives designed for compatibility with all dental materials*, including all BISCO adhesives, and overlaid with a light-cured composite.

Indications For Use

- Base Increment in a Class I or II restoration

Warnings:

- This product may not be appropriate for use in patients who have a history of allergic reactions to methacrylate resins.
- When using dental adhesives, saliva and oral fluid contamination will seriously compromise dentin bonding.

Cautions:

- Cross-contamination: Product may contain items that are designed for one time use. Dispose of used or contaminated accessories. Do not clean, disinfect or reuse.

Precautions:

- Avoid contact with the skin; unpolymerized (meth)acrylate resins may cause skin sensitization in susceptible persons. In case of contact, wash skin with soap and water.
- Working and setting times will be accelerated in the warm intraoral environment.
- To prevent polymerization or clogging, the dual-syringe mixing tip should be left in place until the next application. If used intraorally, wipe mixing tip with disinfectant prior to storage.
- Refer to the manufacturer's instructions. Unlike BISCO Adhesives, not all bonding systems are chemically compatible with self-cured or dual-cured materials.
- Safety data sheet available on request.
- Safety data sheet available at www.Bisco.com.

TECHNICAL INFORMATION

At Room Temperature (20°C/68°F - 25°C/77°F):

Working Time: Minimum 2 min.
Setting Time: Maximum 4 min. 30 sec.

Working time and Setting time may vary based on storage conditions, temperature, humidity, etc.

INSTRUCTIONS FOR USE

1. Isolate tooth and prepare the cavity in a conservative manner. Clean the entire surface with a slurry of pumice and water or pumice with CAVITY CLEANSER™*. Rinse and dry the preparation.
2. Apply an adhesive that is compatible with all dental materials according to manufacturer's instructions.
3. If using individual syringes, dispense by turning syringe handle clockwise. After the desired amount is dispensed, turn handle counterclockwise immediately to prevent excess material extruding from syringe. Mix to a uniform paste (10-15 seconds). To minimize air incorporation, do not over-mix. Utilize instruments designed for composite placement to prevent discoloration of the composite. If using the dual-syringe delivery system, follow the enclosed DUAL-SYRINGE INSTRUCTIONS FOR USE to dispense and deliver **BISFIL 2B**.
4. Using a Low Viscosity Unit-Dose Tip*, place **BISFIL 2B** only to the level of the dentino-enamel junction. Before the end of the working time, place a small increment of a light-cured composite into the preparation. Pack lightly and burnish against the **BISFIL 2B** in order to express entrapped air. Following final set of **BISFIL 2B**, continue to firmly pack a light-cured composite, building incrementally until cavity preparation is filled. Light cure according to manufacturer's instructions.
5. Check the occlusion and begin finishing procedures. After finishing and polishing procedures have been completed, it is recommended to apply a surface sealant, such as FORTIFY™* or BisCover™* LV according to manufacturer's instructions. This will improve marginal integrity and reduce wear and staining.

HYGIENE: Use of commonly available hygienic protective covering to avoid any contamination of the **BISFIL 2B** syringes during treatment is recommended.

DISPOSAL: Refer to community provisions relating to waste. In their absence, refer to national or regional provisions relating to waste.

STORAGE: Store at room temperature (20°C/68°F - 25°C/77°F). Refer to individual component labels for specific expiration dates.

WARRANTY: BISCO, Inc. recognizes its responsibility to replace products if proven to be defective. BISCO, Inc. does not accept liability for any damages or loss, either direct or consequential, stemming from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

¹ Bertolotti R. Posterior Composite Technique Utilizing Directed Polymerization and a Novel Matrix. *Pract Periodont Aesthetic Dent* 1991; 3(4):53-58.

* BISFIL, CAVITY CLEANSER, FORTIFY and BisCover are trademarks of BISCO, Inc.

Low Viscosity Unit-Dose Tips are distributed by BISCO, Inc.

* Refer to manufacturer's instructions. Not all one-bottle adhesives bond to dual-cured or self-cured composites.

BISCO, INC.
1100 W. Irving Park Rd.
Schaumburg, IL 60193
U.S.A.
1-800-247-3368
1-847-534-6000
www.bisco.com



BISICO France
120, allée de la Coudoulette
13680 Lançon De Provence
France
33-4-90-42-92-92