



**39** CLINICAL EVALUATORS

**560** TOTAL USES

**96%** CLINICAL RATING



**Key features:** Resin-modified calcium silicate • Automix • Dual-cured  
• Allows for immediate restoration placement

## Description

**TheraCal PT®** is a biocompatible, dual-cured, resin-modified calcium silicate that is used for pulpotomy treatment.

### TheraCal PT:

- Maintains tooth vitality by performing as a barrier and protectant of the dental pulpal complex.
- Precise and direct placement allows its use in deep cavity preparations.
- The dual-cured set permits immediate placement of the restorative material.
- Alkaline pH of 11.5
- Physical properties that resist breakdown and degradation leading to a durable seal
- Minimum of 45 second working time and maximum of 5-minute set time (dual-cured)
- Secondary indication for direct and indirect pulp capping



## Primary Indication

- Pulpotomies

## Unique Attributes

- Calcium release
- Dual-cured mode allows for a 10-second light cure and immediate placement of restoration
- Moisture tolerant

## Clinical Tip from Evaluators

- Very firm after light curing, so feel confident to place bonding agent and restoration after the 10- second cure.

**“GOOD WORKING TIME AND HARD SET AFTER LIGHT CURING.”**

## Evaluators' Comments

*“Easy application with bendable metal cannula.”*

*“Great consistency - it flows well while adapting to surfaces... even axial walls.”*

*“I believe it is a superior product with good clinical results.”*

*“I already called Bisco and ordered another syringe!”*

*“I needed to apply a lot of pressure to express from the tip.”*

*“Tip was a little wide.”*


## Standout Features

**45+ sec**  
Working Time:  
Minimum 45 seconds at 35°C

**<5 min**  
Setting Time:  
Maximum 5 minutes at 35°C

**1 Visit**  
Time Saver:  
Complete in 1 office visit

## Evaluation Summary:



## Compared to Competitive Products:

