# Bisco

# TransLuma™ Post

## Radiopaque Translucent Fiber Post System

CE 0459

## Instructions for Use

IN-164R5  
Rev. 12/18

TRANSLUMA™™ POST
Radiopaque Translucent Fiber Post System
GENERAL INFORMATION
<p>BISCO's <b>TRANSLUMA POSTS</b> are constructed from unidirectional, prestressed, glass fibers bound in an epoxy resin matrix. This design produces a post with flexural strength exceeding that of metal posts but with a modulus of elasticity similar to that of dentin. The low modulus dissipates stress rather than transferring and concentrating stress in the residual root structures as seen with cast or prefabricated metal posts.</p> <p>Due to its translucency, BISCO's <b>TRANSLUMA POST</b> offers excellent aesthetics for the final restoration without the use of additional masking agents. This post will intimately bond to the tooth structure and core build-up composite while its tapered design respects morphology and conserves tooth structure.</p>
CLINICAL PROCEDURE

With the tooth endodontically treated and filled with gutta-percha, the remaining tooth structure is prepared according to conventional principles

**WARNING:** Drills are sold non-sterile and must be sterilized prior to first use and between subsequent uses in accordance to the instructions provided below.

**CAUTION:** U.S. Federal law restricts this device to sale by or on the order of a dental professional.

- Posts are single use devices and should never be re-used. Structural integrity of the post could be compromised including microscopic damage which could lead to post fracture.

**PRECAUTION:** Gloves must be worn when handling the **TRANSLUMA POSTS** since the post fibers may irritate the skin and the oils on the skin may affect the bond strength of the post. Unpolymerized resins may cause skin sensitization in susceptible persons. Avoid contact with the skin. In case of contact with the skin, wash thoroughly with soap and water. Phosphoric acid/UNI-ETCH is a severe eye irritant and tissue irritant. If accidentally instilled into the eye, flush with copious amounts of water and seek medical attention immediately. In case of contact with other tissues, rinse immediately with plenty of water for several minutes. Injury may result if UNI-ETCH is allowed to remain on the skin or mucosa for extended periods of time.

- Gloves must be worn when handling a post since the post fibers may irritate the skin and the oils on the skin may affect the bond strength of the post.
- The post must be cleaned with alcohol before insertion in the canal.
- Shortening the post should be done outside of the mouth.
- The use of rubber dam is strongly recommended.
- The crown must cover the preparation of at least 1.5 mm of healthy dentin in order to get the ferrule effect.
- Never use pliers to cut the post since the pressure can destroy the structure of the post. In addition, surface alteration of the post (sandblasting or roughening with a bur) is not indicated.
- Safety data sheet available on request.
- Safety data sheet available at [www.bisco.com](#).

**NOTE:** If bonding agents, cementation materials and core materials, other than those identified in these instructions are used, please refer to the manufacturer's instructions for guidance.

- SELECTION OF POST SIZE**

The posts and drills are sized using the ISO Standard numbers that match conventional endodontic files and reamers. Drills are color coded and the packaging for the corresponding posts are marked with the color coding for easy identification as follows:

  - ISO 70 Pre-Shaping Drill: Green
  - ISO 90 Post and Drill: White
  - ISO 100 Post and Drill: Yellow
  - ISO 120 Post and Drill: Blue

- CANAL PREPARATION**
  - General recommendations for post space preparation:
    - The post should be % to % of the length of the root.
    - A minimum of 3-5mm of gutta-percha should remain at the apex.
    - Proper isolation must be used.
  - Remove gutta-percha with a #1 or #2 Peeso Reamer, #3 Gates Glidden or a heated endodontic pluggger.
  - Select the post size according to the ISO number that matches the endodontic files and reamers previously used.
  - Shaping the canal:
    - Use the ISO 70 Pre-Shaping Drill (green) to complete the preliminary preparation. This instrument determines the final depth to which the post will be placed and serves to guide the final preparation with the subsequent drills.
    - An endodontic rubber stopper on the shaft of the drill should be used as a guide for determining the length of the preparation with radiographic verification. To ensure against root perforation while drilling, make sure the drill is always positioned within the canal.
    - Remove debris from the canal with water followed by drying with paper points.
  - Final canal preparation
    - Use the ISO Drills in order, starting with the smallest, until the desired final size is accomplished.

**PLEASE NOTE:** Posts require the sequential use of the drills as follows:

      - ISO 90:** Use the Pre-Shaping Drill ISO 70 (green), followed by the ISO 90 Drill (white).
      - ISO 100:** Use the Pre-Shaping Drill ISO 70 (green), followed by the ISO 90 Drill (white), followed by the ISO 100 Drill (yellow).
      - ISO 120:** Use the Pre-Shaping Drill ISO 70 (green), followed by the ISO 90 Drill (white), followed by the ISO 100 Drill (yellow), then the ISO 120 Drill (blue).
    - Ideally, the final canal preparation allows sufficient space for the post and 30 microns around the post for bonding materials.

- DETERMINING POST LENGTH**
  - Try-in the post and determine the overall length required.
  - Remove the post from the canal and cut with a diamond bur or diamond disc. **WARNING: Never use pliers to cut the post since the pressure can destroy the structure of the post. In addition, surface alteration of the post (sandblasting or roughening with a bur) is contraindicated.**
- POST CEMENTATION**

Please refer to the manufacturer's instructions for guidance when cementing post and building the core.
- FINAL PREPARATION**

Conventional principles of tooth preparation (with proper ferrule) must be followed for optimal results.

  - Preparation for final restoration must end on sound tooth structure.
  - There must be 1.5mm of sound tooth structure beyond the core material.
  - Positive horizontal and vertical walls should exist.
- INFECTION CONTROL**

**Transluma Drills:** these instruments are supplied non-sterile. Disinfect and sterilize instruments before each use. Ultrasonically disinfect instruments with disinfecting or cleaning agents adapted for rotary instruments. Do not disinfect instruments in thermodesinfektor. Inspect instruments visually before each use: change to new drills when they are damaged, corroded or dull. The universal drill may be used only to remove root canal filling partially, and the finishing drills may be used only for post removal canal. Transluma drills must be sterilized individually in an autoclave with the following settings:

  - Wrapped individually in disposable sterilization pouches or tubing that comply with the regulation standard enforced in the country of use (e.g. EN ISO 11607-1).
  - Autoclave: type B complying with regulation standard in force (e.g. EN 13060).
  - Sterilization temperature: 134°C/274°F - sterilization time: 18 minutes.
  - Use a new drill after 12-15 uses.

INFECTION CONTROL FOR TRANSLUMA POST		
	RECOMMENDED	NOT RECOMMENDED
<b>DISINFECTION</b> All chemical disinfection solutions or thermodesinfektor		X
<b>CLEANING</b> Alcohol e.g. ethanol: 1-5 minutes Boiling Ultrasonic cleaning	X	X  X
<b>STERILIZATION</b> (if required) Hot air sterilization Gas sterilization/Chemiclave Steam sterilization/Autoclave (max. 134°C/274°F) 18 minutes, 1 cycle only Wrapped individually in disposable sterilization pouches or tubing that comply with the regulation standard enforced in the country of use (e.g. EN ISO 11607-1). Autoclave: type B complying with regulation standard in force (e.g. EN 13060).	X	X  X

\* Unapproved sterilization may render the post unusable.

**STORAGE:** Store the sterilized components in a dry, dust-free place. If the packaging's integrity appears to be compromised, before using place the component in a new pouch and re-sterilize according to the protocol defined below. The drill should be stored in the sterilization container until required. Containers or pouches must be dry before opening to avoid recontamination of the contents with water. Storage should be in dry, clean conditions at ambient temperature.

**WARNING:** BISCO, Inc. recognizes its responsibility to replace products if proven to be defective. BISCO, Inc. does not accept liability for any damage or loss, either direct or consequential, stemming from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

\* TRANSLUMA is a trademark of BISCO, Inc.

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