Instructions for Use

BISCO, Inc.
1100 W. Irving Park Road
Schaumburg, IL 60193
U.S.A.
1-847-534-6000
1-800-247-3368

IN-029R9
Rev. 10/17

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed Dentist.

RESINOMER™
Amalgam Bonding/Luting System

GENERAL INFORMATION

RESINOMER is a dual-cured composite with viscosity and flow characteristics that make it ideal for bonded amalgam restorations, as a dental liner and as a luting cement. The hydrophilic monomers used in RESINOMER are based on ALL-BOND® 2 chemistry and promote adhesion to dental metals. Chemical adhesion occurs between the metallic surface and diisobutylate dimethylacrylate (DSDM), one of the monomers present in RESINOMER. The metal bonding capability provides an exceptional environment for bonded amalgam restorations. A bonding medium bond is created as amalgam is condensed into the uncured RESINOMER layer. In addition to adhesive bonding, RESINOMER provides a highly effective seal against microleakage. RESINOMER is intended for use with adhesives designed for compatibility with all dental materials, including BISCO adhesives. RESINOMER's bond to ALL-BOND 2 produces resin adhesion that exceeds the cohesive strength of dentin.

Indications for Use:

- 1. Bonded Amalgam Restorations
- 2. Liner/Base under Direct/Indirect Restorations
- 3. Cementation of Metal-Based Restorations
- 4. Cementation of Fiber and Metallic Posts
- 5. Bonding Orthodontic Brackets

Warnings:

- This product may not be appropriate for use in patients who have a history of allergic reactions to methacrylic resins.
- When using dental adhesives, saliva and oral fluid contamination will seriously compromise dentin bonding.
- Cautions:**
- Cross-contamination: Product may contain items that are designed for one time use. Dispose of used or contaminated accessories. Do not clean, disinfect or reuse.
- Precautions:**
- Avoid contact with the skin; unpolymerized resin may cause skin sensitization in susceptible persons. In case of contact, wash skin with soap and water.
- Store in refrigerator when not in use. Allow refrigerated RESINOMER to reach room temperature before use.
- Work and setting times will be accelerated in the warm intraoral environment.
- If using ALL-BOND 2, cement will set quickly once it comes in contact with the primers on the tooth.
- To prevent polymerization or clogging, the dual-syringe lipp should be left in place until the next application. If left intraorally, wipe excess lip with disinfector prior to storage.
- See individual component labels for specific expiration dates.
- Use BISCO adhesives in bonding systems that are chemically compatible with self-cured or dual-cured materials. Refer to the adhesive manufacturer's instructions.
- Do not use hydrogen peroxide, EDTA, or eugenol containing materials prior to the adhesive procedure as these materials may inhibit bonding.
- Safety data sheet available on request.
- Safety data sheet available at www.Bisco.com.

TECHNICAL INFORMATION

Individual Syringes:

Working Time: Minimum 3 min. (At Room Temperature, 20°C/68°F - 25°C/77°F)

Setting Time: Maximum 3 min. 20 sec. (Intraoral, 37°C/98.6°F)

Dual-Syringes:

Working Time: Minimum 2 min. (At Room Temperature, 20°C/68°F - 25°C/77°F)

Setting Time: Maximum 3 min. 20 sec. (Intraoral, 37°C/98.6°F)

Working time and Setting time may vary based on storage conditions, temperature, humidity, etc.

NOTE: Due to the unique chemistry of RESINOMER, refrigeration is necessary. Allow refrigerated RESINOMER to reach room temperature before use.

NOTE: RESINOMER left on the mixing pad may appear not to set. This is due to oxygen inhibition of the polymerization process. Rest assured, RESINOMER will polymerize (set) in a normal fashion. The user can confirm this by placing a mylar strip over the excess material on the pad; it will set within the specified period of time.

INSTRUCTIONS FOR USE

BONDED AMALGAM RESTORATIONS

1. Use a matrix band with non-stick coating, or rub wax over the area that will be exposed to RESINOMER to prevent adhesion to the band.

2. Isolate tooth and prepare the cavity in a conservative manner. Clean the entire surface with a slurry of powder and a cleaner such as CAVITY CLEANSER™ or water.

3. Etch using an etchant, such as UNI-ETCH® w/BAC or ETCH-37™ w/BAC, according to manufacturer's instructions.

4. Apply an adhesive that is compatible with all dental materials according to manufacturer's instructions.

5. If using the dual-syringe delivery system, follow the enclosed DUAL-SYRINGE INSTRUCTIONS FOR USE for dispensation and delivery of RESINOMER. If using individual syringes, mix equal amounts of RESINOMER and catalyst into a uniform paste (10-15 seconds) until the paste is a uniform color.

6. Light cure RESINOMER to the floor or desired internal area of the preparation.

7. Trifurture and condense amalgam into RESINOMER-inlaid preparation using traditional condensation procedures. Do not be concerned with displacement of RESINOMER.

8. Slightly overlap, dent and burnish amalgam over margins. Begin denting and after curing to enhance the bond and seal the margins. A slight flash of amalgam onto the enamel surface is desirable and may be unavoidable due to the adhesive qualities of RESINOMER.

9. OPTIONAL: RESINOMER is double-toothable. To expedite the polymerization, light cure for 20 seconds.

NOTE: RESINOMER left on the mixing pad may appear not to set. This is due to oxygen inhibition of the polymerization process. Rest assured, RESINOMER will polymerize (set) in a normal fashion. The user can confirm this by placing a mylar strip over the excess material on the pad; it will set within the specified period of time.

LAYER/BASE UNDER DIRECT/INDIRECT RESTORATIONS

1. Isolate tooth and prepare the cavity in a conservative manner. Clean the entire surface with a slurry of powder and a cleaner such as CAVITY CLEANSER™ or water.

2. Etch using an etchant, such as UNI-ETCH® w/BAC or ETCH-37™ w/BAC, according to manufacturer's instructions.

3. Apply an adhesive that is compatible with all dental materials according to manufacturer's instructions.

4. If using the dual-syringe delivery system, follow the enclosed DUAL-SYRINGE INSTRUCTIONS FOR USE for dispensation and delivery of RESINOMER. If using individual syringes, mix equal amounts of RESINOMER base and catalyst into a uniform paste (10-15 seconds) until the paste is a uniform color.

5. Apply a thin coat of RESINOMER to the floor or desired internal area of the preparation.

6. Remove excess with a brush or instrument. Complete clean-up during the initial set or gel state.

FIBER AND METALLIC POST CEMENTATION

A. Post Space Preparation

1. Etch using a etchant, such as UNI-ETCH w/BAC or ETCH-37 w/BAC, according to manufacturer's instructions.

2. Etch using an etchant, such as UNI-ETCH w/BAC or ETCH-37 w/BAC, according to manufacturer's instructions.

3. Apply an adhesive that is compatible with all dental materials according to manufacturer's instructions.

4. Sandblast the internal surfaces of the restoration according to manufacturers instructions.

5. Seal the restoration with passive pressure.

6. Remove excess with a brush or instrument. Complete clean-up during the initial set or gel state.

B. Cementation

1. Sandblast the internal surfaces of the restoration, rinse and dry.

2. Apply a metal primer to the internal surface of the restoration according to manufacturers instructions.

3. Using the dual-syringe delivery system, follow the enclosed DUAL-SYRINGE INSTRUCTIONS FOR USE for dispensation and delivery of RESINOMER. If using individual syringes, mix equal amounts of RESINOMER base and catalyst into a uniform paste (10-15 seconds) until the paste is a uniform color.

4. Fill the internal surface of the restoration with RESINOMER.

5. Seal the restoration with passive pressure.

6. Apply an adhesive that is compatible with all dental materials according to manufacturer's instructions.

7. Proceed with restorative treatment.

METAL-BASED RESTORATION CEMENTATION

A. Tooth Preparation

1. Remove provisional restorations; clean the preparation and try-in the definitive restoration.

2. Etch using an etchant, such as UNI-ETCH w/BAC or ETCH-37 w/BAC, according to the manufacturer's instructions.

3. Apply an adhesive that is compatible with all dental materials according to manufacturer's instructions.

4. Sandblast the internal surfaces of the restoration, rinse and dry.

5. Seal the restoration with passive pressure.

6. Remove excess with a brush or instrument. Complete clean-up during the initial set or gel state.

C. Scellement

1. Sandblast the internal surfaces of the restoration, rinse and dry.

2. Apply a metal primer to the internal surface of the restoration according to manufacturers instructions.

3. Using the dual-syringe delivery system, follow the enclosed DUAL-SYRINGE INSTRUCTIONS FOR USE for dispensation and delivery of RESINOMER. If using individual syringes, mix equal amounts of RESINOMER base and catalyst into a uniform paste (10-15 seconds) until the paste is a uniform color.

4. Fill the canal dry with a paper point; repeat with a new paper point until the paper point returns dry from the canal. This step is important for removal of any pooled primers to prevent interference in complete sealing of the post.

5. Light cure as necessary, according to manufacturer's instructions.

6. Coat the post with adhesive and light cure, according to manufacturer's instructions.

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B. Post Cementation

1. Coat the post with adhesive and light cure, according to manufacturer's instructions.

2. If using the dual-syringe delivery system, follow the enclosed DUAL-SYRINGE INSTRUCTIONS FOR USE for dispensation and delivery of RESINOMER. If using individual syringes, mix equal amounts of RESINOMER base and catalyst into a uniform paste (10-15 seconds) until the paste is a uniform color.

3. Using the dual-syringe delivery system, follow the enclosed DUAL-SYRINGE INSTRUCTIONS FOR USE for dispensation and delivery of RESINOMER. If using individual syringes, mix equal amounts of RESINOMER base and catalyst into a uniform paste (10-15 seconds) until the paste is a uniform color.

4. Coat the apical end of the post with adhesive and light cure, according to manufacturer's instructions.

5. Light cure as necessary, according to manufacturer's instructions.

6. Coat the apical end of the post with adhesive and light cure, according to manufacturer's instructions.

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