**Bisco**

CE 0459

# ORTHO Bracket Paste LC

## Light-Cured Bracket Bonding System

## Instructions for Use

BISCO, Inc.  
1100 W. Irving Park Rd.  
Schaumburg, IL 60193  
U.S.A.  
1-847-534-6000  
1-800-247-3368

IN-088R6  
Rev. 10/17

**CAUTION:** U.S. Federal law restricts this device to sale by or on the order of a dental professional.

### ORTHO BRACKET PASTE LC\* Light-Cured Bracket Bonding System

#### GENERAL INFORMATION

**ORTHO BRACKET PASTE LC** is a light-cured orthodontic direct bonding system. This system combines a micromechanical bond to enamel and a chemical bond between ONE-STEP+ and **ORTHO BRACKET PASTE LC** to generate maximum bonding strength from all types of dentition to ceramic, metal or plastic brackets.

#### Indications For Use

Bracket Bonding

#### Warnings:

- The phosphoric acid in the etchant is a severe eye and skin irritant. Injury may result if the etchant is allowed to remain on the skin or mucosa for extended periods of time. If accidentally splashed into the eye, flush with copious amounts of water and seek medical attention immediately. In case of contact with other tissues, rinse immediately with plenty of water for several minutes.
- When using dental adhesives, saliva contamination will seriously compromise dentin bonding.
- ONE-STEP is highly flammable.

#### Cautions:

- Cross-contamination: Product may contain items that are designed for one time use. Dispose of used or contaminated unit-dose components and accessories. Do not clean, disinfect or reuse.
- Do not use prophylactic paste which may contain oils, since it may leave a film that will reduce etching quality and decrease bonding strength.
- Do not rub the syringe tip over the enamel to avoid fracture of the enamel rods.

#### Precautions:

- Keep etchant away from direct sunlight.
- Avoid contact with the skin; unpolymerized resins may cause skin sensitization in susceptible persons. In case of contact, wash skin with soap and water.
- See individual component labels for specific expiration dates.
- Safety data sheet available on request.
- Safety data sheet available at [www.Bisco.com](http://www.Bisco.com).

#### INSTRUCTIONS FOR USE

- Prophylaxis - Thoroughly clean the teeth to be bonded with pumice in a rubber cup or bristle brush. Carefully avoid irritating the gingival tissue since slight bleeding may cause contamination and could affect bonding quality. **CAUTION:** Do not use prophylactic paste which may contain oils, since it may leave a film that will reduce etching quality and decrease bonding strength.
- Rinse thoroughly with plenty of water. Isolate the teeth and dry with the air syringe.
- Apply ETCH-37™ directly on the enamel surface, limited to the area to be bonded for 30 seconds (up to 60 seconds for deciduous teeth). **CAUTION:** Do not rub the syringe tip over the enamel to avoid fracture of the enamel rods. Alternatively, use an absorbent pellet to carefully dab the acid onto each tooth to achieve proper acid penetration. Also do not etch for a period over 90 seconds, which may cause the formation of insoluble calcium phosphate salt crystals which may decrease bonding strength.
- Rinse thoroughly with plenty of water for about 5 seconds on each tooth. Dry the teeth with oil-free air. The etched areas should appear dull and frosty white. If not, re-etch areas as needed.
- Hold bottle at a 45° angle and dispense 1-2 drops of ONE-STEP into a mixing well. Using a brush, apply a thin coat of ONE-STEP to the etched enamel.
- Using a gentle stream, gently but thoroughly air dry to evaporate solvent.
- Apply **ORTHO BRACKET PASTE LC** to the back of the bracket. Place the bracket on the sealed tooth, gently pressing it to its final position. Remove the excess material. Be careful to not displace all the paste from underneath the bracket. With the light probe as close as possible to the bracket, light cure the bracket for 30 seconds, following the directions below:
  - Metal brackets** - light cure for 20 seconds from the incisal and 10 seconds from the gingival or mesial surfaces.
  - Ceramic or Plastic brackets** - light cure for 30 seconds through the archwire slot from the facial surface.
- Once the paste has been properly cured, an active archwire can be placed immediately.

**BRACKET REMOVAL:** **ORTHO BRACKET PASTE LC** can be easily removed with pin and ligature cutting pliers or other de-bonding devices. Tooth surfaces may be further cleaned with scaling instruments and polished with rotary instruments.

**HYGIENE:** Re-cap syringes with the luer lock cap. Use of commonly available hygienic protective covering to avoid any contamination of syringes during treatment is recommended.

**DISPOSAL:** Refer to community provisions relating to waste. In their absence, refer to national or regional provisions relating to waste.

**STORAGE:** Store at room temperature (20°C/68°F - 25°C/77°F). Refer to individual component labels for specific expiration dates.

**WARRANTY:** BISCO, Inc. recognizes its responsibility to replace products if proven to be defective. BISCO, Inc. does not accept liability for any damages or loss, direct or consequential, stemming from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

\* ONE-STEP is a registered trademark of BISCO, Inc.

ORTHO BRACKET PASTE LC is manufactured by BISCO, Inc.

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ETCH-37 is a trademark of BISCO, Inc.
## BISCO, INC.1100 W. Irving Park Road Schaumburg, IL 60193 U.S.A.