**BISCO**

**CE 0459**

**AELITE™ Composite** *Light-Cured*

**Instructions for Use**

**BISCO, Inc.**  
 1100 W. Irving Park Road  
 Schaumburg, IL 60193  
 U.S.A.  
 1-800-247-3368  
 1-847-534-6000  
 www.bisco.com

**BISCO France**  
 208, allée de la Coudoulette  
 13680 Lançon de Provence  
 France  
 33-4-90-42-92-92

IN-131R12  
 Rev. 9/21

**Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed Dentist.**

**en** **ÆLITE™ Composites**

**GENERAL INFORMATION**  
 The **ÆLITE** system of composites consists of light-cured resin-based dental restorative materials that provide exceptional aesthetics for all classes of direct restorations. This system is designed to accommodate the "single-composite restoration" technique and the "multi-layering" technique.

**ÆLITE Aesthetic Enamel - Reinforced Nanofil Composite**  
**Aesthetic Enamel** is a light-cured, fluorescing, highly reinforced, nanofil composite. Due to its reinforced composition, **Aesthetic Enamel** can be used for all classes of direct restorations as the enamel layer. The nano-sized filler particles contribute to **Aesthetic Enamel's** excellent polishability while the reinforcing filler allows for high filler loading giving **Aesthetic Enamel** maximum strength and durability.

**ÆLITE All-Purpose Body - Universal Microhybrid Composite**  
**All-Purpose Body** is a light-cured, radiopaque, universal microhybrid composite. **All-Purpose Body** is designed to blend with surrounding dentition when used alone or in conjunction with **ÆLITE Aesthetic Enamel**.

**ÆLITE LS Posterior - Low Shrinkage (LS) Posterior Composite**  
**ÆLITE LS Posterior** is a light-cured, highly-filled, radiopaque, low shrinkage posterior hybrid composite. Due to the high filler content, **LS Posterior** aids in establishing occlusal contours and interproximal contacts. In the layering technique, **Aesthetic Enamel** may be placed over the **LS Posterior**.

- Indications for use:**  
 Principle uses of **ÆLITE Aesthetic Enamel**:
- Enamel layer for Class I, II, III, IV, V, restorations
  - Distal-most veneers
  - Direct veneers
  - Composite and porcelain repairs
  - Splinting periodontally involved teeth

- Principle uses of **ÆLITE All-Purpose Body**:
- Class I, II, III, IV, V, restorations
  - Direct aesthetic veneers for missing stains and anomalies
  - Core build-ups to replace missing tooth structure
  - Composite and porcelain repairs
  - Lingual increment of anterior restorations
  - To mask up discolored adjacent teeth at the base of a preparation (opaque shades)

- Principle uses of **ÆLITE LS Posterior**:
- Class I, II and VI restorations
  - Core build-ups to replace missing tooth structure

- Warnings:**
- If accidentally splashed into the eye, flush with copious amounts of water and seek medical attention immediately. In case of contact with other tissues, rinse immediately with plenty of water for several minutes.

- Cautions:**
- Product may contain items that are designed for one time use. Dispose of used or contaminated tools. Do not clean, disinfect or reuse.
  - Use of commonly available hygienic protective covering to avoid any contamination of syringes during treatment is recommended. If syringe becomes contaminated, discard. Do not clean or disinfect.

- Precautions:**
- Avoid contact with the skin; unpolished methacrylate resins may cause skin sensitization in susceptible persons. In case of contact, wash skin with soap and water.
  - When using dental adhesives, saliva contamination will seriously compromise dentin bonding.
  - Refer to individual component labels for specific expiration dates.
  - Safety data sheet available on request.
  - Safety data sheet available at www.bisco.com.

INSTRUCTIONS FOR USE	
1. Select the desired type and shade of <b>ÆLITE</b> composite prior to isolating and preparing the tooth: <b>LS Posterior</b> : for Class I, Class II and Class VI <b>All-Purpose Body</b> : Universal Composite for Class I - Class VI <b>Aesthetic Enamel</b> : Enamel Layer for Class I - Class VI	
<b>Tip:</b> To aid in shade selection, place a small amount of material on the tooth and cure. The composite is easily removed with the use of an instrument since it has not been bonded to the tooth.	
2. Clean the preparation with a slurry of pumice and water. Rinse and dry the preparation.	
3. Isolate the preparation with the <b>ÆLITE</b> composite. Apply a 1 mm layer of <b>TheraCal LC™</b> directly on and around the exposed pulp, then light cure for 20 seconds. If prep is near the pulp, apply a 1 mm layer of <b>TheraCal LC™</b> as liner and light cure for 20 seconds.	
4. Place matrix, and apply an adhesive according to the manufacturer's instructions.	
5. Dispense desired amount of <b>ÆLITE</b> composite onto a light-protected mixing pad by turning the syringe handle clockwise. (Tip: To prevent excessive material from extruding, turn handle counterclockwise immediately after the desired amount is dispensed.) If using Uni-Dose™ Tips, inject <b>ÆLITE</b> composite directly into preparation.	
<b>Precautions:</b>	
<ul style="list-style-type: none"> <li>Uni-dose tips are for single use only. Do not clean, disinfect or reuse.</li> </ul>	
<b>NOTE:</b> The properties of composite naturally tend itself to adhere (stick) to your instrument. To overcome this tackiness, <b>MODELING RESIN™</b> was developed to aid in the placement and shaping of your composites. If lightly wetting your instrument with <b>MODELING RESIN</b> , any of your direct composite restorations can be easily placed and manipulated (sculpted).	
6. Using the selected shades, place 1-2 mm increments of composite into the cavity preparation. Light cure each increment with a curing light unit (refer to curing chart for curing time intervals). Keep the light as close to the composite as possible without contact. Continue to build incrementally until the cavity preparation is filled to the occlusal margin and the restoration top is close to the desired shade. Each surface of the final increment should be light-cured for the amount of time specified in the curing chart.	

Product	For each 1-2 mm increment, light-cure for:	For the final increment, light-cure each surface for:
AESTHETIC ENAMEL	20 seconds	30 seconds
ALL-PURPOSE BODY Opaque shades<1mm	20 seconds 30 seconds	30 seconds
LS POSTERIOR A1, A2, A3, B1, C2, D3	20 seconds	40 seconds
LS POSTERIOR A3.S, C4	40 seconds	40 seconds

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7. Check the occlusion and finish. To minimize or eliminate polishing, apply **BiCover™ LV** Liquid Polish surface sealants such as **FORTIFY™** or **FORTIFY PLUS** may also be used according to manufacturer's instructions. The use of any of these products will seal the margins and reduce composite wear and staining.

**HYGIENE:** Use of commonly available hygienic protective covering to avoid any contamination of syringes during treatment is recommended.  
**DISPOSAL:** Refer to community provisions relating to waste. In their absence, refer to national or regional provisions relating to waste.  
**WARNING:** **BISCO, Inc.** recognizes its responsibility to replace products proven to be defective. **BISCO, Inc.** does not accept liability for any damages or loss, either direct or consequential, stemming from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

- TheraCal** is a registered trademark of **BISCO, Inc.**
- MODELING RESIN** is manufactured by **BISCO, Inc.**
- ÆLITE, BiCover** and **FORTIFY** are trademarks of **BISCO, Inc.**

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• TherCa/L er registert varemärke som tilhører BISCO, Inc.
• MODELING RESIN produseres av BISCO, Inc.
• TherCa/L er een geadponeerd handelsmerk van BISCO, Inc.
• MODELING RESIN wordt vervaardigd door BISCO, Inc.
• ELITE, BioCover en FORTIFY zijn handelsmerken van BISCO, Inc.