BISCO

TheraCal PT®Dual-Cured

Dual-Cured Resin-Modified Calcium Silicate Pulpotomy Treatment

Instructions for Use

 BISCO, Inc.
1100 W. Irving Park Road
Schaumburg, IL 60193
U.S.A.
1-800-247-3368
1-847-534-6000
www.bisco.com

IN-201-ZH
Rev. 8/26

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed Dentist.

en	Theracal PT**
GENERAL INFORMATION	
TheraCal PT is a biocompatible, dual-cured, resin-modified calcium silicate that is used to treat exposed dentin. TheraCal PT's precise placement allows its use in deep cavity preparations. The dual-cured set permits immediate placement of the restorative material. Its proprietary formulation allows for a command set with a light curing unit while maintaining ease of placement due to thixotropic properties.	
Model specifications	Description
H-34112P-zh	Including TheraCal PT: H-34010 1Syringe(4g) AUTO-MIX CANNULA TIPS: X-81269P 15Mixing Tips
H-34113P-zh	Including TheraCal PT: H-34111 1Syringe(2g) AUTO-MIX CANNULA TIPS: X-81271P 8Mixing Tips

Specification:
Name of registrant/manufacturer: BISCO, Inc.
Residence/production address: 1100 W. Irving Park Road Schaumburg, IL 60193 U.S.A.
Contact details: 1-800-247-3368 / 1-847-534-6000 / www.bisco.com
Name of agent/after-sales service unit:Waterside Medical Group
Address of agent/after-sales service unit:East Half of 2nd Floor and 5th Floor, Building 3, Zone 3, No.20 Datong Road, High-Tech Zone, Suzhou, Jiangsu, China
Tel: 13818791340
Medical Device Registration No.:20263170262
Product Technical Requirement No.:20263170262
Revision date of the instructions:2026.7
Product performance, main structure composition
This product is resin-modified calcium silicate.

Scope of application / Indications
1.Pulpotomy
2.Direct pulp capping
3.Indirect pulp capping and cavity lining.

Contraindication
TheraCal PT is not recommended for use with patients who have a history of severe allergic reaction to methacrylate resins.

Warnings
• TheraCal PT contains polymerizable methacrylate monomers. Avoid prolonged or repeated exposure with skin, oral soft tissues, and eyes. Irritation and possible corneal damage may result. Skin rash, oral mucosa irritation, or other allergic reactions (allergic contact dermatitis) may result in susceptible individuals. Do not take internally.
• TheraCal PT is not recommended for use with patients who have a history of severe allergic reaction to meth acrylate resins.

Cautions
Cross-contamination:
• Product may contain items that are designed for one time use. Dispose of used or contaminated tips. Do not clean, disinfect or reuse.
• Use of commonly available hygienic protective covering to avoid any contamination of syringes during treatment is recommended. If syringe becomes contaminated, discard. Do not clean or disinfect.

Precautions
• The success of pulp capping procedures may be influenced by preparation contamination. Use of a rubber dam or adequate isolation is recommended.
• Patient history of pain influences appropriateness of usage.
• Pulp hemorrhage must be controlled prior to the application of TheraCal PT. If hemorrhage cannot be controlled, consideration should be given for the initiation of endodontic therapy over direct pulp capping procedures.
• Protective eyewear should be worn by the patient and dental staff.
• TheraCal PT is a dual-cured material. Proceed immediately once material has been dispensed.
• To prevent polymerization or clogging, the dual-syringe mixing tip should be left in place until the next application. If used intraorally, wipe mixing tip with disinfectant prior to storage.
• When chlorhexidine and sodium hypochlorite are used together staining can occur. If used together there should be a saline rinse between.
• Procedures that de-vitalize the pulp (for example, formocresol) may prevent secondary bridge formation.
• See individual component labels for specific expiration dates.
• Safety data sheet available on request.
• Safety data sheet available at www.bisco.com.
• Instructions for Use available upon request.
• Instructions for Use available at www.bisco.com.

Technical Information
Working Time: Minimum 45 seconds at 35°C
Setting Time: Maximum 5 minutes at 35°C
Working time and setting time may vary based on storage conditions (i.e. temperature, humidity, etc.).

INSTRUCTION FOR USE

TheraCal PT Mixing Procedures
1) Remove the syringe cap from the dual-syringe or mixing tip.
2) Dispense a small amount of material onto a mixing pad to eliminate any voids in each chamber of the dual-syringe.
3) Attach a mixing tip to the dual-syringe by lining up the key and keyway. Then turn the locking mechanism clockwise.
4) Pressing the plunger will mix and dispense the TheraCal PT. (See Technical Information section for working and setting time information.)

Pulpotomy
1) Under rubber dam isolation, remove the infected dentin with a round bur and/or a hand excavator.
2) Gain access to the pulp chamber and remove the infected coronal pulp using sterile burs or instruments.
3) If there is bleeding in the pulp, first achieve hemostasis by placing a sterile cotton pellet dampened with sterile saline (or preferred alternative solution).
4) Adapt a matrix around the tooth if a proximal wall is missing.
5) Place TheraCal PT directly in the pulp chamber, and ensure good adaptation to the cavity walls and margins.
6) Light cure TheraCal PT for 10 seconds.
7) Place desired adhesive, base, and/or restoration following manufacturer's directions for restoring tooth.

Pulp Exposures (Direct Pulp Capping)
1) Under rubber dam isolation, complete cavity preparation.
2) First achieve hemostasis by placing a sterile cotton pellet dampened with sterile saline (or preferred alternative solution) gently onto the exposure site.
3) Gently blot preparation with sterile cotton pellet. Leave visibly moist.
4) Apply TheraCal PT directly to the exposed pulp. Cover all the exposed areas and extend TheraCal PT at least 1mm onto sound dentin surrounding the exposure.
5) Light cure TheraCal PT for 10 seconds.
6) Place desired adhesive, base, and/or restoration following manufacturer's directions. Continue restoring tooth.

Protective Liner and Base for use under a variety of substrates
1) Under rubber dam isolation, perform conventional cavity preparation. Remove all infected carious tooth structure. Leave preparation visibly moist.
2) Apply the TheraCal PT directly to the cavity floor. Manipulate into a smooth surface covering all deep dentin areas.
3) Light cure TheraCal PT for 10 seconds.

Hygiene: Use of commonly available hygienic protective covering to avoid any contamination of the TheraCal PT syringes during treatment is recommended.

Storage: Store in a dry place at room temperature (20°C/68°F - 25°C/77°F). TheraCal PT syringe should remain in the pouch until ready to use. Once pouch is opened, use TheraCal PT within 90 days.
Transport environment: Regular shipping, avoid direct sunlight, rain and snow.

Date of manufacture: See the label

Validity period: See the label

* BISCO and TheraCal PT are trademarks of BISCO, Inc., registered in the U.S. Patent and Trademark Office.

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