



All-Bond Universal®
Light-Cured Dental Adhesive

Universal Primer™
Dual-Cured Adhesive

All-Bond SE®
Self-Etching Adhesive

All-Bond 3®
Universal Dental Adhesive System

All-Bond 2®
Universal Dental Adhesive

**One-Step®/
One-Step® Plus**
Universal Dental Adhesive/
Filled Universal Dental Adhesive

BISCO Adhesives

Light-Cured



Dual-Cured



Compatible with
Light-, Dual-, and
Self-Cured Materials



Total-Etch Technique



Self-Etch Technique



Selective-Etch Technique



No Additional Activator
Required for Compatibility
with Dual-Cured Products



Single Bottle



Compatible with
Both Direct & Indirect
Restorations



Storage at
Room Temperature



MDP Containing



Best Uses

A light-cured, single bottle with no additional activator for all direct and indirect restorations to be used with any bonding technique.

A dual-cure primer/adhesive designed for those who prefer to not light-cure the adhesive layer under indirect restorations.

A light-cured adhesive for all direct and indirect restorations to be used in the self-etch bonding technique.

A dual-cured, hydrophobic adhesive for all direct and indirect restorations to be used in the total-etch bonding technique.

A dual-cure, hydrophilic, BPDM containing monomer adhesive system to be used in the total-etch bonding technique.

A light-cured, single bottle with no additional activator for all direct and indirect restorations to be used in the total-etch bonding technique.

Ideal for those who want ONE adhesive in the office for all restorations.

*Use One-Step Plus for a "Filled" adhesive option to reduce post-operative sensitivity.

