

5 REASONS TO TRY

Quantium™

Universal Composite

1 Durable Restorations

Quantium's high strength, wear resistance, and low shrinkage create durable restorations.

2 Low Shrinkage

Quantium features low volumetric shrinkage, which minimizes microleakage and ensures marginal adaptation.

3 Excellent Esthetics

Quantium is designed to be easily polished and stain resistant to create and preserve esthetic results and maintain their appearance over time.

4 Easy to Handle

Quantium delivers smooth, non-sticky handling and its unit dose capsules can be warmed for an even greater flowable experience, depending on the clinician's needs and preference.

5 Highly Radiopaque

Quantium is highly radiopaque for easy identification on radiographs aiding in post-treatment evaluation, allowing clinicians to confidently monitor the integrity of the restoration.



We are here to help! Give us a call: 1-800-247-3368 or visit us online: www.bisco.com


BISCO
Rx Only MC-50525QU